

## **Product Usage Guide**

## **All Purpose Sanitiser**


Health & Safety Please refer to relevant MSDS

An all-purpose sanitiser based on natural fermentation extracts and mild, biodegradable surfactant. Kills bacteria to EN1276. Approved to EN14476 for enveloped viruses, including coronavirus 2019 n-CoV / SARS-CoV-2.










Tip bottle and use guide to measure correct dilution in the measuring chamber. Remove cap and pour.




Dilute 10ml of product into a 750ml trigger spray. 1L makes up to 100 Spray Bottles.



Spray surface, leave for the required contact time, (15 seconds for enveloped viruses, 30 seconds for bacteria) and wipe.\*\*





Tip bottle and use guide to measure correct dilution in the measuring chamber. Remove cap and pour.



Dilute 60ml of product into a 5L bucket.

Please note: Colour of product is for demonstration purposes. Red line shows correct dosing level.



Mop surface, leave for the required contact time, (15 seconds for enveloped viruses, 30 seconds for bacteria). No rinsing necessary. Leave to air dry.\*\*





\*\*please ensure surfaces are clean prior to sanitising

