

Product Usage Guide

All Purpose Sanitiser

Health & Safety Please refer to relevant MSDS An all-purpose sanitiser based on natural fermentation extracts and mild, biodegradable surfactant. Kills bacteria to EN1276. Approved to EN14476 for enveloped viruses, including coronavirus 2019 n-CoV / SARS-CoV-2.









Dilute 8ml of product into a 750ml trigger spray. 5L makes up to 625 spray bottles.






Dilute 60ml of product into a 5L bucket.



Spray surface, leave for the required contact time, (15 seconds for enveloped viruses, 30 seconds for bacteria) and wipe.**



Mop surface, leave for the required contact time, (15 seconds for enveloped viruses, 30 seconds for bacteria). No rinsing necessary. Leave to air dry.**



**please ensure surfaces are clean prior to sanitising

Rinsing not required after application



Pelican Pump Sizes 30ml

For further information and MSDS contact: BioHygiene, Pantglas Industrial Estate, Bedwas, Caerphilly, CF83 8DR T: 029 2067 4094 E: info@biohygiene.co.uk www.biohygiene.co.uk